

Drosophila Vials

I. APPLICATIONS

1. Storage accessory for Drosophila related research

II. FEATURES

1. Materials: polystyrene
2. Glass-like clarity for easy sample viewing
3. Disposable

III. VIAL DESIGN

1. Available in Narrow (25 x 95mm) and Wide (28.5 x 95mm)
2. Both narrow and wide vials are available in bulk and tray packaging

IV. PACKAGING

1. Bulk, 500 pieces/case
2. Tray, 100 vials/tray, 500 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
51-0500	Narrow, 25 x 95mm vial	Bulk, 500 pieces/case	No
51-0501	Narrow, 25 x 95mm vial	Tray, 500 pieces/case	No
51-0800	Wide, 28.5 x 95mm vial	Bulk, 500 pieces/case	No
51-0801	Wide, 28.5 x 95mm vial	Tray, 500 pieces/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C
51-0500	Drosophila Vial	95	O.D. 22.6	O.D. 25
51-0800	Drosophila Vial	95	O.D. 25.4	O.D. 28.5



51-0500



51-0800