

Drosophila Bottle

I. APPLICATIONS

1. Storage accessory for Drosophila related research

II. FEATURES

1. Materials: polypropylene
2. Autoclavable
3. Disposable
4. Transparent for easy sample viewing

III. BOTTLE DESIGN

1. 6oz (177ml) bottle
2. Molded graduations in 25ml increments up to 100ml
3. Square bottom
4. Available in both bulk and tray packaging

IV. PACKAGING

1. Bulk, 200 pieces/case
2. Tray, 25 bottles/tray, 200 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
51-17720	Drosophila bottles, square bottom	Bulk, 200 pieces/case	No
51-17721	Drosophila bottles, square bottom	Tray, 200 pieces/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H
51-17720	Drosophila Bottle	103.3	74.6	12.0	6.0	O.D. 43.0	O.D.38.5	55.6	34.2



51-17720