

Drosophila Plugs

I. APPLICATIONS

1. Fit for drosophila vials and bottles for drosophila related research

II. FEATURES

1. Materials: high density foam
2. Disposable
3. Offer minimal shedding as well as reduce evaporation
4. Superior barrier to mite infestation
5. Stretches yet maintains tight fit after repeated removal and closing
6. Flat surface accepts marking pen

III. PLUG DESIGN

1. Available in Narrow and Wide
2. Exact fit with narrow 0.9 in. (25mm) or wide 1.1 in. (28.5mm) diameter vials or 6oz (177ml) bottles

IV. PACKAGING

1. Bulk, 100 pieces/bag, 10 bags/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
51-0505	Fit for 25 x 95mm vial	Bulk, 10,000 pieces/case	No
51-0805	Fit for 28.5 x 95mm vial	Bulk, 10,000 pieces/case	No
51-17725	Fit for 60oz bottle	Bulk, 10,000 pieces/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


IX. PRODUCT PARAMETER

Dimensions=mm


Cat #	Description	A	B
51-0505	Drosophila Plug	28	45
51-0805	Drosophila Plug	34	45
51-17725	Drosophila Plug	45	45



51-0505



51-0805



51-17725